

Specification For Xylene

TEST	LIMITS
Assay (by GC) (Xylene Isomers & Ethyl Benzene)	98.5% min
Paraxylene	Per Lot
Metaxylene	Per Lot
Orthoxylene	Per Lot
Ethyl Benzene	25% max.
Boiling Range	137-143°C
Color (APHA)	10 max.
Appearance	Clear
Residue After Evaporation	10 ppm max.
Substances Darkened By H_2SO_4	To Pass
Sulfur Compounds	1 ppm max.
Water	0.05% max.
Specific Gravity @ 15.56°C	0.865-0.875
Benzene Content	0.01% max.
Toluene Content	0.5% max.
Total C9 Aromatics	Per Lot
Non-Aromatics	0.30% max.

Packing: - 200 Litres refilled in used Drums

Disclaimer: The information contained herein with regards to all products specification, is based on technical data and tests which we believe to be reliable. However, we assume no responsibility for results obtained or damage, if any, incurred through the application of data given, since conditions of use are outside our control.